# PUSH™

InBody

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Quick Start Guide **USA**

Quick Start Guide

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---

# PUSH™

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# **PUSH**<sup>TM</sup>

## **Quick Start Guide** **USA**

- This product uses ultrasonic waves to measure height easily and quickly.
- Please read the 'Quick Start Guide' before using the product.

## 1. Product Components

The PUSH consists of the following components.  
Please make sure all of the components are present.

### ❶ PUSH



### ❷ Battery – AAA(1.5V) X4


### ❸ Quick Start Guide

## 2. Appearance and Part Names

### 3. Battery Safety Instructions

Insert 4 AAA batteries so the positive (+) and negative (-) terminals align correctly.

- 1) Insert batteries so the positive (+) and negative (-) terminals align correctly.  
Batteries may rupture or leak if inserted incorrectly.
- 2) Do not mix old and new batteries. Only use AAA batteries from the same brand.  
Replace all 4 AAA batteries from your device at the same time.
- 3) Only use alkaline or manganese batteries.
- 4) Do not handle batteries with wet hands.
- 5) Remove batteries from devices that aren't being used.
- 6) Dispose of batteries according to local regulations.
- 7) Batteries may damage your device or cause a fire if used improperly.

## 4. InBody hi App

You can monitor and manage your child's growth conveniently with InBody hi App.

Go to the Google Play Store or Apple App Store and download the InBody hi App.

< Android OS >

< iOS >


## 5-1. Measuring Height with 'PUSH'

- 1) Turn on the PUSH by pressing the button. Indoor temperature and humidity level will be checked automatically.



- 2) Press the button again to change to measuring mode. Ultrasonic waves will measure the distance to the ground. Since ultrasonic waves spread vertically & horizontally, make sure there are no obstacles like walls and/or toys nearby.

※ Make sure the person being measured keeps his/her arms by his/her sides. This will help make sure the ultrasonic waves are not being obstructed.

- 3) Place the end of the PUSH that doesn't have the ultrasonic wave sensor on the top of the head. Balance the PUSH vertically and horizontally until "PUSH" is displayed on the screen.

- 4) Press the button to take measurement. You can re-measure the height by pressing the button when you see horizontal bar on the screen.

"Please press the button 

- ※ To turn it off, press the button for 3 seconds or let it idle for 20 seconds.

## 5-2. Working with App to Measure Height

- 1) Launch the InBody hi App.
- 2) Select the 'hi' button located in the center of the bottom navigation on the home screen.
- 3) Select the user you want to record in the App and press the Next button.
- 4) Select 'InBody PUSH' and press the Next button.
- 5) Press the button on your PUSH device to measure the height.
  - \* Please refer to '5-1 Measuring height with PUSH' for how to measure height.
  - \* If the height measured with PUSH does not display in the App, please check if your smartphone's Bluetooth is turned on.
- 6) Select the measured height value.
  - \*Up to 10 height values will be displayed on the App screen.
- 7) Enter the weight to complete.
- 8) (Optional) After entering the weight, you can upload photos and record a Growth diary in the Growth album.
- 9) You can check the measured results in hi App > Growth Development.

## 6. Error Messages

Following messages will appear on the LCD screen if:

PUSH is not balanced horizontally.

PUSH will show up on the screen when the PUSH is balanced.

**E050** The person being measured is shorter than 1ft. 7.7in. PUSH cannot measure people shorter than 1ft. 7.7in.

**E999** Ultrasonic waves were not detected. Please clear the obstacles, including carpet and rug. For best results, use the PUSH on a hard surface like wood or tile.

**E200** The person being measured is taller than 6ft. 6.7in. PUSH cannot measure people taller than 6ft. 6.7in.

**E000** The sensor is covered. Please make sure the ultrasonic wave sensor is not covered by objects like head or hands.

**bAtt** Battery icon will appear on the screen when the battery is low. 'bAtt' will appear on the screen when the battery is too low to measure. The device will turn off automatically. Replace batteries.

## 7. Cautions

- 1) Keep the device out of the reach of children. They might get hurt if they throw or drop the device.
- 2) Do not drop or put too much pressure on the ultrasonic wave sensor.
- 3) Take measurements with the PUSH in similar environments for consistent results.
- 4) Do not measure height on carpets or blankets. Height will not be measured accurately.
- 5) Take measurements away from obstacles like tables, toys and walls for accurate results.
- 6) Keep the device away from water.
- 7) If the PUSH is arbitrarily disassembled, or damaged by natural disaster, InBody will void the warranty.

## 8. Specifications

Measurement Item	Height, Temperature, Humidity	Height Range	1ft 7.7 in ~ 6ft 6.7 in (50 ~ 200cm)
Error Range	± 0.2 in (± 0.5cm)	Measurement Duration	1 sec. (Max 7 sec.)
Sensor	Ultrasound, Temperature, Humidity	Other Measuring Range	Temperature 60.8 ~ 104°F (± 33.8°F) 16 ~ 40°C (± 1°C)
Current Consumption	Below 100mA		
Power Consumption	DC 3V (1.5V AAA Battery 4EA)		Humidity 30 ~ 80% (± 7%)
Dimension	1.69(W) X 1.73(L) X 11.8(H) : inch	Display Type	LCD
	43(W) x 44(L) 300(H) : mm	Device Weight	0.52 lb (238g)
Operation Environment	60.8 ~ 104°F (16 ~ 40°C), 30 ~ 80%RH	Storage Environment	14 ~ 158°F (-10 ~ 70°C)

\* Specifications are subject to change without notice.

Only representatives from InBody service center are allowed to repair and inspect the PUSH. Please contact the Product Support team at [Support@inbody.com](mailto:Support@inbody.com).

## 9. Product Warranty

### Product Warranty Standards

- Warranty service will be effective only within the warranty period provided and only for an inherent manufacturer's defect.
- Any provisions not specified herein will follow the Fair Trade Commission notice of consumer dispute resolution standards.
- Only trained InBody personnel are allowed to repair and inspect PUSH.

### Warranty period

- One year from the date of purchase
- If the purchase date cannot be verified, the warranty period will become effective three months after manufacture date.

Types of Consumer Damages		Within the warranty period	After expiry of the warranty period
If major malfunction occurs within the first 30 days of purchase due to defect.		Exchange or Repair free of charge	Repair at a charge
Repairs Under Warranty	Any malfunction occurring due to a defect within 1 year of the purchase date	Repair or exchange with refurbished product of equal or better condition	

## **Paid Service**

Any of the following cases will void the product warranty.

- Arbitrary disassembly of the product.
- Unauthorized modification of product controls.
- Repair by a person other than the repair representative designated by InBody.
- Failure caused by the consumer's fault:
  - Product failure and damage caused by the user's negligence (fall, moisture intrusion, shock, damage, forcible control, etc.).
  - Product failure and damage caused by the violation of the Quick Start Guide.
- Others:
  - Product failure caused by a natural disaster (e.g., fire, flood, abnormal power, etc.).
  - Worn-out or used consumables (batteries or others).
- Shipping costs for repair or replacement may be required at the cost of the consumer

# **PUSH<sup>TM</sup>**

## **Quick Start Guide**

- This product uses ultrasonic waves to measure height easily and quickly.
- Please read the 'Quick Start Guide' before using the product.

# 1. Product Components

The PUSH consists of the following components.  
Please make sure all of the components are present.

## ❶ PUSH



## ❷ Battery – AAA(1.5V) X4


## ❸ Quick Start Guide

# 2. Appearance and Part Names

### 3. Battery Safety Instructions

Insert 4 AAA batteries so the positive (+) and negative (-) terminals align correctly.

- 1) Insert batteries so the positive (+) and negative (-) terminals align correctly.  
Batteries may rupture or leak if inserted incorrectly.
- 2) Do not mix old and new batteries. Only use AAA batteries from the same brand.  
Replace all 4 AAA batteries from your device at the same time.
- 3) Only use alkaline or manganese batteries.
- 4) Do not handle batteries with wet hands.
- 5) Remove batteries from devices that aren't being used.
- 6) Dispose of batteries according to local regulations.
- 7) Batteries may damage your device or cause a fire if used improperly.

## 4. InBody hi App

You can monitor and manage your child's growth conveniently with InBody hi App.

Go to the Google Play Store or Apple App Store and download the InBody hi App.

< Android OS >

< iOS >


## 5-1. Measuring Height with 'PUSH'

- 1) Turn on the PUSH by pressing the button. Indoor temperature and humidity level will be checked automatically.



- 2) Press the button again to change to measuring mode. Ultrasonic waves will measure the distance to the ground. Since ultrasonic waves spread vertically & horizontally, make sure there are no obstacles like walls and/or toys nearby.

※ Make sure the person being measured keeps his/her arms by his/her sides. This will help make sure the ultrasonic waves are not being obstructed.

- 3) Place the end of the PUSH that doesn't have the ultrasonic wave sensor on the top of the head. Balance the PUSH vertically and horizontally until "PUSH" is displayed on the screen.

- 4) Press the button to take measurement. You can re-measure the height by pressing the button when you see horizontal bar on the screen.

"Please press the button 

- ※ To turn it off, press the button for 3 seconds or let it idle for 20 seconds.

## 5-2. Working with App to Measure Height

- 1) Launch the InBody hi App.
- 2) Select the 'hi' button located in the center of the bottom navigation on the home screen.
- 3) Select the user you want to record in the App and press the Next button.
- 4) Select 'InBody PUSH' and press the Next button.
- 5) Press the button on your PUSH device to measure the height.
  - \* Please refer to '5-1 Measuring height with PUSH' for how to measure height.
  - \* If the height measured with PUSH does not display in the App, please check if your smartphone's Bluetooth is turned on.
- 6) Select the measured height value.
  - \* Up to 10 height values will be displayed on the App screen.
- 7) Enter the weight to complete.
- 8) (Optional) After entering the weight, you can upload photos and record a Growth diary in the Growth album.
- 9) You can check the measured results in hi App > Growth Development.

## 6. Error Messages

Following messages will appear on the LCD screen if:

PUSH is not balanced horizontally.

**PUSH** will show up on the screen when the PUSH is balanced.

**E050** The person being measured is shorter than 50cm. PUSH cannot measure people shorter than 50cm.

**E999** Ultrasonic waves were not detected. Please clear the obstacles, including carpet and rug. For best results, use the PUSH on a hard surface like wood or tile.

**E200** The person being measured is taller than 200cm. PUSH cannot measure people taller than 200cm.

**E000** The sensor is covered. Please make sure the ultrasonic wave sensor is not covered by objects like head or hands.

**bAtt** Battery icon will appear on the screen when the battery is low. 'bAtt' will appear on the screen when the battery is too low to measure. The device will turn off automatically. Replace batteries.

## 7. Cautions

- 1) Keep the device out of the reach of children. They might get hurt if they throw or drop the device.
- 2) Do not drop or put too much pressure on the ultrasonic wave sensor.
- 3) Take measurements with the PUSH in similar environments for consistent results.
- 4) Do not measure height on carpets or blankets. Height will not be measured accurately.
- 5) Take measurements away from obstacles like tables, toys and walls for accurate results.
- 6) Keep the device away from water.
- 7) If the PUSH is arbitrarily disassembled, or damaged by natural disaster, InBody will void the warranty.

## 8. Specifications

Measurement Item	Height, Temperature, Humidity	Height Range	50 ~ 200cm (1ft 7.7 in ~ 6ft 6.7 in)
Error Range	± 0.5cm (± 0.2 in)	Measurement Duration	1 sec. (Max 7 sec.)
Sensor	Ultrasound, Temperature, Humidity	Other Measuring Range	Temperature 16 ~ 40°C (± 1°C) 60.8 ~ 104°F (± 33.8°F)
Current Consumption	Below 100mA		
Power Consumption	DC 3V (1.5V AAA Battery 4EA)		Humidity 30 ~ 80% (± 7%)
Dimension	43(W) x 44(L) 300(H) : mm	Display Type	LCD
	1.69(W) X 1.73(L) X 11.8(H) : inch	Device Weight	238g (0.52 lb)
Operation Environment	16 ~ 40°C (60.8 ~ 104°F), 30 ~ 80%RH	Storage Environment	-10 ~ 70°C (14 ~ 158°F)

\* Specifications are subject to change without notice.

Only representatives from InBody service center are allowed to repair and inspect the PUSH.

## 9. Product Warranty

### Product Warranty Standards

- Warranty service will be effective only within the warranty period provided and only for an inherent manufacturer's defect.
- Any provisions not specified herein will follow the Fair Trade Commission notice of consumer dispute resolution standards.
- Only trained InBody personnel are allowed to repair and inspect PUSH.

### Warranty period

- One year from the date of purchase
- If the purchase date cannot be verified, the warranty period will become effective three months after manufacture date.

Types of Consumer Damages		Within the warranty period	After expiry of the warranty period
If major malfunction occurs within the first 30 days of purchase due to defect.		Exchange or Repair free of charge	Repair at a charge
Repairs Under Warranty	Any malfunction occurring due to a defect within 1 year of the purchase date	Repair or exchange with refurbished product of equal or better condition	

## **Paid Service**

Any of the following cases will void the product warranty.

- Arbitrary disassembly of the product.
- Unauthorized modification of product controls.
- Repair by a person other than the repair representative designated by InBody.
- Failure caused by the consumer's fault:
  - Product failure and damage caused by the user's negligence (fall, moisture intrusion, shock, damage, forcible control, etc.).
  - Product failure and damage caused by the violation of the Quick Start Guide.
- Others:
  - Product failure caused by a natural disaster (e.g., fire, flood, abnormal power, etc.).
  - Worn-out or used consumables (batteries or others).
- Shipping costs for repair or replacement may be required at the cost of the consumer

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TEL: +31-20-238-6080 FAX: +31-6-5734-1858

Website: nl.inbody.com E-mail: info.eu@inbody.com

Website (CS): csinbody.eu E-mail(CS): support@inbody.com

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Website: de.inbody.com E-mail: erfolg@inbody.com

**InBody UK [UNITED KINGDOM]**

Orega, Belmont Road, Uxbridge, UB8 1HE, United Kingdom

TEL: +44-1530-569620

Website: [uk.inbody.com](http://uk.inbody.com) E-mail: [uk@inbody.com](mailto:uk@inbody.com)

**InBody Oceania [AUSTRALIA]**

Main office: Level 8, 1 York Street, SYDNEY, NSW 2000, Australia

Showroom: U2/82-86 Minnie Street, Southport, Queensland

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**InBody Asia [MALAYSIA & SINGAPORE]**

Unit 3A-11, Oval Damansara, 685 Jalan Damansara Kuala Lumpur,  
WP KL 60000 Malaysia

TEL: +60-3-7732-0790

FAX: +60-3-7733-0790

Website: [inbodyasia.com](http://inbodyasia.com)

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E-mail: [info.mx@inbody.com](mailto:info.mx@inbody.com)

**InBody India [INDIA]**

57/57 A, 1st Floor, Raj Industrial Complex, Military Road, Marol,  
Andheri (East). Mumbai- 400059, Maharashtra, India

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This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions : (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operations.

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